

Product Datasheet

Mouse APOB ELISA Kit (orb438012)

Description	Mouse Apolipoprotein B (APOB) Sandwich ELISA Kit
Reactivity	Mouse
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	E9Q414
Sensitivity	1.22ng/mL
Entrez	238055
Expiration Date	Please enquire.