

Product Datasheet

Rat APOB ELISA Kit (orb438014)

Description	Rat Apolipoprotein B (APOB) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q7TMA5
Sensitivity	11.3ng/mL
Entrez	54225
Expiration Date	Please enquire.