

Product Datasheet

Human APOC3 ELISA Kit (orb438021)

Description	Human Apolipoprotein C3 (APOC3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P02656
Sensitivity	0.138ng/mL
Entrez	345
Expiration Date	Please enquire.