

Product Datasheet

Human APOC4 ELISA Kit (orb438024)

Description	Human Apolipoprotein C4 (APOC4) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P55056
Sensitivity	37pg/mL
Entrez	346
Expiration Date	Please enquire.