

Product Datasheet

Mouse APOC4 ELISA Kit (orb438025)

Description	Mouse Apolipoprotein C4 (APOC4) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q61268
Sensitivity	0.17ng/mL
Entrez	11425
Expiration Date	Please enquire.