

Product Datasheet

Human APOD ELISA Kit (orb438026)

Description	Human Apolipoprotein D (APOD) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P05090
Sensitivity	0.64ng/mL
Entrez	347
Expiration Date	Please enquire.