

Product Datasheet

Bovine APOH ELISA Kit (orb438034)

Description	Bovine Apolipoprotein H (APOH) Sandwich ELISA Kit
Reactivity	Bovine
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P17690
Sensitivity	1.87ng/mL
Entrez	281006
Expiration Date	Please enquire.