

Product Datasheet

Mouse APOH ELISA Kit (orb438036)

Description	Mouse Apolipoprotein H (APOH) Sandwich ELISA Kit
Reactivity	Mouse
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, saliva and other biological fluids.
Sample Types	serum, plasma, saliva and other biological fluids.
Uniprot ID	Q01339
Sensitivity	1.42ng/mL
Entrez	11818
Expiration Date	Please enquire.