

Product Datasheet

Human APP ELISA Kit (orb438045)

Description	Human Amyloid Precursor Protein (APP) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cerebrospinal fluid and other biological fluids.
Sample Types	serum, plasma, cerebrospinal fluid and other biological fluids.
Uniprot ID	P05067
Sensitivity	35pg/mL
Entrez	351
Expiration Date	Please enquire.