

Product Datasheet

Rat APP ELISA Kit (orb438047)

Description	Rat Amyloid Precursor Protein (APP) Sandwich ELISA Kit
Reactivity	Rat
Range	1.25-80ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P08592
Sensitivity	0.49ng/mL
Entrez	54226
Expiration Date	Please enquire.