

Product Datasheet

Human AQP2 ELISA Kit (orb438050)

Description	Human Aquaporin 2, Collecting Duct (AQP2) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, urine and other biological fluids.
Sample Types	tissue homogenates, urine and other biological fluids.
Uniprot ID	P41181
Sensitivity	0.122ng/mL
Entrez	359
Expiration Date	Please enquire.