

Product Datasheet

Mouse ARC ELISA Kit (orb438061)

Description	Mouse Activity Regulated Cytoskeleton Associated Protein (ARC) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9WV31
Sensitivity	0.058ng/mL
Entrez	11858
Expiration Date	Please enquire.