

Product Datasheet

Bovine AREG ELISA Kit (orb438063)

Description	Bovine Amphiregulin (AREG) Sandwich ELISA Kit
Reactivity	Bovine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	A5PJE7
Sensitivity	12.8pg/mL
Entrez	538751
Expiration Date	Please enquire.