

Product Datasheet

Rat AREG ELISA Kit (orb438066)

Description	Rat Amphiregulin (AREG) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P24338
Sensitivity	12.7pg/mL
Entrez	29183
Expiration Date	Please enquire.