

Product Datasheet

Rat ARHGEF7 ELISA Kit (orb438073)

Description	Rat Rho Guanine Nucleotide Exchange Factor 7 (ARHGEF7) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O55043
Sensitivity	0.064ng/mL
Entrez	114559
Expiration Date	Please enquire.