

Product Datasheet

Rat ARNTL ELISA Kit (orb438078)

Description	Rat Aryl Hydrocarbon Receptor Nuclear Translocator Like Protein (ARNTL) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9EPW1
Sensitivity	0.057ng/mL
Entrez	29657
Expiration Date	Please enquire.