

Product Datasheet

Mouse ARNT ELISA Kit (orb438079)

Description	Mouse Aryl Hydrocarbon Receptor Nuclear Translocator (ARNT) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P53762
Sensitivity	0.054ng/mL
Entrez	11863
Expiration Date	Please enquire.