

Product Datasheet

Human ARPC2 ELISA Kit (orb438082)

Description	Human Actin Related Protein 2/3 Complex Subunit 2 (ARPC2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O15144
Sensitivity	0.068ng/mL
Entrez	10109
Expiration Date	Please enquire.