

Product Datasheet

Mouse ARPC4 ELISA Kit (orb438083)

Description	Mouse Actin Related Protein 2/3 Complex Subunit 4 (ARPC4) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P59999
Sensitivity	0.043ng/mL
Entrez	68089
Expiration Date	Please enquire.