



Product Datasheet

Blood Group ABH antibody (orb43809)

Description	Mouse Monoclonal to Blood Group ABH.
Reactivity	Human
Conjugation	Other
Tested Applications	IHC-P
Immunogen	Mixture of erythrocytes of group A1 and glycoprotein fraction isolated from saliva of secretors with blood group A.
Target	Blood Group ABH
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (paraffin sections): The antibody HE-10 is excellent as a tumour marker in patients of blood group B and O. Flow cytometry: Recommended dilution: 1-3 : 100.
Isotype	Mouse IgM
Clonality	Monoclonal
Purity	Hybridoma culture supernatant concentrated by membrane ultrafiltration.
Dilution Range	Immunohistochemistry (paraffin sections): The antibody HE-10 is excellent as a tumour marker in patients of blood group B and O. Flow cytometry: Recommended dilution: 1-3 : 100.
Expiration Date	12 months from date of receipt.



Separation of human erythrocytes from bl...



Flow cytometry surface staining pattern ...