

Product Datasheet

Human ASTL ELISA Kit (orb438109)

Description	Human Astacin Like Metallo Endopeptidase (ASTL) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q6HA08
Sensitivity	0.060ng/mL
Entrez	431705
Expiration Date	Please enquire.