

Product Datasheet

Blood Group A antibody (orb43811)

Description	Mouse Monoclonal to Blood Group A.
Reactivity	Human
Conjugation	Other
Tested Applications	AGG
Immunogen	Mixture of erythrocytes of group A1 and glycoprotein fraction isolated from saliva of secretors with blood group A.
Target	Blood Group A
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgM
Clonality	Monoclonal
Purity	Hybridoma culture supernatant concentrated by membrane ultrafiltration.
Expiration Date	12 months from date of receipt.

Clone	Agglutination of tested blood groups				
	A:	B	A:B	A:B	O
HE-183 (anti-A)	+	-	+	+	-
HE-185 (anti-A)	+	-	+	+	-
HEB-29 (anti-B)	-	+	+	+	-
HE-24 (anti-A:B)	+	-	+	-	-

Agglutination of particular blood groups...