

Product Datasheet

Human ATG16L1 ELISA Kit (orb438115)

Description	Human Autophagy Related Protein 16 Like Protein 1 (ATG16L1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q676U5
Sensitivity	0.063ng/mL
Entrez	55054
Expiration Date	Please enquire.