

Product Datasheet

Mouse ATM ELISA Kit (orb438120)

Description	Mouse Ataxia Telangiectasia Mutated (ATM) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q62388
Sensitivity	5.6pg/mL
Entrez	11920
Expiration Date	Please enquire.