

Product Datasheet

Mouse ATP1b3 ELISA Kit (orb438125)

Description	Mouse ATPase, Na ⁺ /K ⁺ Transporting Beta 3 Polypeptide (ATP1b3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P97370
Sensitivity	0.066ng/mL
Entrez	11933
Expiration Date	Please enquire.