

Product Datasheet

Rat ATP1b3 ELISA Kit (orb438126)

Description	Rat ATPase, Na ⁺ /K ⁺ Transporting Beta 3 Polypeptide (ATP1b3) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q63377
Sensitivity	0.071ng/mL
Entrez	25390
Expiration Date	Please enquire.