

Product Datasheet

Human ATP1b4 ELISA Kit (orb438127)

Description	Human ATPase, Na ⁺ /K ⁺ Transporting Beta 4 Polypeptide (ATP1b4) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UN42
Sensitivity	0.106ng/mL
Entrez	23439
Expiration Date	Please enquire.