

## Product Datasheet

### Rat ATP1b4 ELISA Kit (orb438128)

|                          |                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Rat ATPase, Na <sup>+</sup> /K <sup>+</sup> Transporting Beta 4 Polypeptide (ATP1b4) Sandwich ELISA Kit |
| <b>Reactivity</b>        | Rat                                                                                                     |
| <b>Range</b>             | 0.156-10ng/mL                                                                                           |
| <b>Form/Appearance</b>   | Traditional                                                                                             |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt.                |
| <b>Note</b>              | For research use only                                                                                   |
| <b>Application notes</b> | tissue homogenates and other biological fluids.                                                         |
| <b>Sample Types</b>      | tissue homogenates and other biological fluids.                                                         |
| <b>Uniprot ID</b>        | <b>Q9R193</b>                                                                                           |
| <b>Sensitivity</b>       | 0.055ng/mL                                                                                              |
| <b>Entrez</b>            | <b>84396</b>                                                                                            |
| <b>Expiration Date</b>   | Please enquire.                                                                                         |