

Product Datasheet

Human ATP2A2 ELISA Kit (orb438129)

Description	Human ATPase, Ca++ Transporting, Cardiac Muscle, Slow Twitch 2 (ATP2A2) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P16615
Sensitivity	0.126ng/mL
Entrez	488
Expiration Date	Please enquire.