

Product Datasheet

Mouse ATP2A2 ELISA Kit (orb438130)

Description	Mouse ATPase, Ca++ Transporting, Cardiac Muscle, Slow Twitch 2 (ATP2A2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O55143
Sensitivity	0.059ng/mL
Entrez	11938
Expiration Date	Please enquire.