

Product Datasheet

Mouse ATP4b ELISA Kit (orb438132)

Description	Mouse ATPase, H ⁺ /K ⁺ Exchanging Beta Polypeptide (ATP4b) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P50992
Sensitivity	0.066ng/mL
Entrez	11945
Expiration Date	Please enquire.