

Product Datasheet

Rat ATP6AP2 ELISA Kit (orb438133)

Description	Rat ATPase, H ⁺ Transporting, Lysosomal Accessory Protein 2 (ATP6AP2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	075787
Sensitivity	0.122ng/mL
Entrez	10159
Expiration Date	Please enquire.