

Product Datasheet

Human ATP7b ELISA Kit (orb438134)

Description	Human ATPase, Cu++ Transporting Beta Polypeptide (ATP7b) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P35670
Sensitivity	0.056ng/mL
Entrez	540
Expiration Date	Please enquire.