

Product Datasheet

Human aZGP1 ELISA Kit (orb438141)

Description	Human Alpha-2-Glycoprotein 1, Zinc Binding (aZGP1) Sandwich ELISA Kit
Reactivity	Human
Range	4.687-300ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, urine, saliva, seminal plasma, sweat and other biological fluids.
Sample Types	serum, plasma, urine, saliva, seminal plasma, sweat and other biological fluids.
Uniprot ID	P25311
Sensitivity	1.9ng/mL
Entrez	563
Expiration Date	Please enquire.