

Product Datasheet

Human bACE1 ELISA Kit (orb438150)

Description	Human Beta-Site APP Cleaving Enzyme 1 (bACE1) Sandwich ELISA Kit
Reactivity	Human
Range	0.25-16ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cerebrospinal fluid and other biological fluids.
Sample Types	serum, plasma, cerebrospinal fluid and other biological fluids.
Uniprot ID	P56817
Sensitivity	0.11ng/mL
Entrez	23621
Expiration Date	Please enquire.