

Product Datasheet

Mouse bACE1 ELISA Kit (orb438151)

Description	Mouse Beta-Site APP Cleaving Enzyme 1 (bACE1) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cerebrospinal fluid and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cerebrospinal fluid and other biological fluids.
Uniprot ID	P56818
Sensitivity	5.9pg/mL
Entrez	23821
Expiration Date	Please enquire.