

Product Datasheet

Human BAFF ELISA Kit (orb438153)

Description	Human B-Cell Activating Factor (BAFF) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9Y275
Sensitivity	0.055ng/mL
Entrez	10673
Expiration Date	Please enquire.