

Product Datasheet

Human BAFFR ELISA Kit (orb438156)

Description	Human B-Cell Activation Factor Receptor (BAFFR) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q96RJ3
Sensitivity	0.061ng/mL
Entrez	115650
Expiration Date	Please enquire.