

Product Datasheet

Human BANP ELISA Kit (orb438161)

Description	Human BTG3 Associated Nuclear Protein (BANP) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8N9N5
Sensitivity	0.055ng/mL
Entrez	54971
Expiration Date	Please enquire.