

Product Datasheet

Bovine Bax ELISA Kit (orb438164)

Description	Bovine Bcl2 Associated X Protein (Bax) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	002703
Sensitivity	0.124ng/mL
Entrez	280730
Expiration Date	Please enquire.