

Product Datasheet

Human BCAR1 ELISA Kit (orb438170)

Description	Human Breast Cancer Anti-Estrogen Resistance 1 (BCAR1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P56945
Sensitivity	0.053ng/mL
Entrez	9564
Expiration Date	Please enquire.