

Product Datasheet

Rat BCHE ELISA Kit (orb438172)

Description	Rat Butyrylcholinesterase (BCHE) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9JKC1
Sensitivity	0.44ng/mL
Entrez	65036
Expiration Date	Please enquire.