

Product Datasheet

Rat BCKDK ELISA Kit (orb438173)

Description	Rat Branched Chain Alpha-Ketoacid Dehydrogenase Kinase (BCKDK) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q00972
Sensitivity	0.043ng/mL
Entrez	29603
Expiration Date	Please enquire.