

Product Datasheet

Human BCL2L11 ELISA Kit (orb438174)

Description	Human Bcl2 Like Protein 11 (BCL2L11) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O43521
Sensitivity	0.114ng/mL
Entrez	10018
Expiration Date	Please enquire.