

Product Datasheet

Mouse BFP ELISA Kit (orb438186)

Description	Mouse Brain Finger Protein (BFP) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q96DY5
Sensitivity	0.066ng/mL
Entrez	22671
Expiration Date	Please enquire.