

Product Datasheet

Human BHLHB8 ELISA Kit (orb438189)

Description	Human Basic Helix Loop Helix Domain Containing Protein B8 (BHLHB8) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q7RTS1
Sensitivity	0.071ng/mL
Entrez	168620
Expiration Date	Please enquire.