

Product Datasheet

Human BIRC7 ELISA Kit (orb438196)

Description	Human Baculoviral IAP Repeat Containing Protein 7 (BIRC7) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q96CA5
Sensitivity	0.073ng/mL
Entrez	79444
Expiration Date	Please enquire.