

Product Datasheet

Rat BLVRA ELISA Kit (orb438201)

Description	Rat Biliverdin Reductase A (BLVRA) Sandwich ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P46844
Sensitivity	0.231ng/mL
Entrez	116599
Expiration Date	Please enquire.