

Product Datasheet

Rat BMP7 ELISA Kit (orb438216)

Description	Rat Bone Morphogenetic Protein 7 (BMP7) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	G3V6W8
Sensitivity	6.1pg/mL
Entrez	85272
Expiration Date	Please enquire.