

Product Datasheet

Mouse BMPER ELISA Kit (orb438220)

Description	Mouse BMP Binding Endothelial Regulator (BMPER) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8CJ69
Sensitivity	0.054ng/mL
Entrez	73230
Expiration Date	Please enquire.